



ANALYZE PRECISELY
MARKET SCIENTIFICALLY



WOOD'S



UV



RGB



PPL



CPL

SKIN ANALYZER MC88

THE AESTHETICIAN PREFERRED DEVICES



CE RoHS

PRODUCT INTRODUCTION Features

1



Unique Design



Clear 15 Pictures



Easy to Operate



Free Training



Visual Analysis



Marketing Products

The MEICET MC88 Skin Image analyzer is a software and hardware integrated system that utilizes image analysis and processing technology.

It is designed to assist in observing skin texture, pigmentation, and skin barrier. The system features five spectral photography modes, including RGB light, cross-polarized light, parallel-polarized light, UV light, and Wood's light. Based on these five spectra, the system captures five corresponding spectral images.

The system analyzes these five spectral images using algorithmic techniques to generate a total of 15 images. These images, along with the final analysis report, assist beauty professionals in conducting comprehensive and accurate analysis of facial skin conditions.



Golden



Black

◉ Show with iPad



◉ Side display



◉ UV light display



◉ Hidden



*Apple iPad is not included. We do not have the right to sell iPad. Thanks for your understanding.



Spectra	RGB Light, Cross-polarized Light,UV Light, Parallel-polarized Light, Wood's Light
Lighting Technology	Solid-State LED
Adapter	Input: AC100~240V, 50-60Hz, 1.2A Output: DC 24V, 2A
Power Consumption	40 W
Net Weight	10 kg
Gross Weight	17 kg
Main Unit Dimensions(mm)	400 (L)*430 (W)*550 (H)
Wireless technology and Working Frequency	BLE, 2402MHZ-2480MHZ
Maximum transmission power of wireless function	4dbm

*iPad Specification Requirements for Selection

iPad System Requirements: iOS System
Storage Memory Requirements: Recommended 128GB and above



What demands can

MEICET MC88 HELP YOU FULFILL?

Core Advantage



Clear 15 Pictures

Provide clear and accurate skin images



Marketing Products


Increase sales revenue



Visual Analysis

Reveal Hidden Skin Problems



Promote your store's products or services: such as moisturizing, sun protection, spot removal, anti-aging, whitening, cleansing, and repair, to effectively increase your sales.



These 15 images output by MC88 system can help to detect different skin problems, such as sensitive skin, enlarged pores, uneven skin tone, pigmentation, wrinkles, porphyrin, skin texture, inflammation, etc.

Simultaneous Comparison of Skin Symptoms ASSISTANCE WITH ANALYTICAL FEATURES


Comparison



Compare different skin symptom images of the same time, to find out truth of skin problems.

Comparison of Identical Skin Symptoms at Different Times BEFORE-AFTER COMPARISON

Comparison




Sensitivity-Heatmap




Sensitivity-Red Area



Spots Removing



Wrinkles




Grid function + Comparison function

Compare the same skin symptom images of different time, to present products effect and get customers' trust, With help of grid function, the effect of tightening and lifting can be checked.


Increase the exposure of the store and products MARKETING YOUR PRODUCTS

Marketing

7




Recommended products on the report




These reports can be printed out or sent to customers' email directly so that the exposure of your store and products can be increased, and the impression of customers can be deepened, thereby growing store visibility and product sales.

Tagging




By directly annotating skin issues on the image, effective visual analysis can be conducted.



View images and reports across multiple devices. Go beyond geographical limitations

Customize Report

"FREE LOGO REPLACEMENT" AND "HOME PAGE CAROUSEL IMAGES IN THE APP"



Free logo replacement



Setting Banners

When exporting reports, you can customize the logo according to your needs. Additionally, on the app, you can replace the promotional banner based on your recent requirements.

Customize Report WATERMARK SETTINGS


10

☒ Time watermark

Watermark types

Time watermark

Watermark position



☒ Text watermark

input "Your text"

Watermark types

Text watermark

Watermark content

MECET

Watermark position



☒ Original image exporting


input "Null"

Watermark types

Text watermark

Watermark content

Watermark position




Added watermark feature with three setting options: Time watermark, Text watermark, and Original image exporting. Effectively enhances brand impression and strengthens copyright protection.

Additionally, it is possible to set the watermark position, effectively avoiding important detection areas.


synchronized with the current device lock screen password

DATA SPACE ENCRYPTION

① "Enable data space encryption"



② "View customer's historical testing information"



③ "Enter the current device lock screen password (sync)"





SHANGHAI MAY SKIN INFORMATION TECHNOLOGY CO., LTD



6th Floor, Building 7, No. 7001, Zhongchun Road, Minhang District,
Shanghai 201100, China



www.meicet.com



info@meicet.com



+86 131 6722 3337